

# ACR Smartphone Urinalysis

Home urine testing kit & application



## How does it work?

The Albumin to Creatinine Ratio test is built around the existing ACR point of care dipstick. The test is suitable for lay users and can be performed anywhere using the ACR test kit and application.

ACR testing is used for the management and early detection of chronic kidney disease.

Results are available immediately and can easily be shared with a clinician or integrated into an existing electronic medical record.

## Performance

- Semi-quantitative analysis of micro albumin (10 - 150 mg/L), creatinine (10 - 300 mg/dL) and albumin to creatinine ratio (mg/g).
- Results are returned following the KDIGO albumin:creatinine thresholds (A1, A2, A3)
- Full adherence service provided including delivery to patient's home, proven to significantly increase adherence to recommended albuminuria screenings.
- CE Marked for sale in the European Union.

## Content

- 1 x Albumin:Creatinine dipstick
- 1 x urine receptacle and lid
- 1 x Colour-Board
- User manual
- Smartphone application can easily be downloaded from the app store

