

Dip.io Smartphone Urinalysis

Home urine testing kit & application



How does it work?

The Dip test uses a standard urine dipstick, and serves as an indication for a wide range of conditions including urinary tract infections, ketosis, health in pregnancy as well as diseases of the kidney and liver.

The test is suitable for lay users and can be performed anywhere using the Dip kit and application.

Results are available immediately and can easily be shared with a clinician or integrated into an existing electronic medical record.

Performance

- Semi-Quantitative analysis of 10 parameter dipstick indicating: Leukocytes, Nitrites, Glucose, Ketones, Protein, Blood, pH, Urobilinogen, Bilirubin, Specific Gravity
- Analytical sensitivity equivalent to point of care semi-quantitative urinalysis devices.
- CE Marked for sale in the European Union and FDA class II 510K clearance.
- Application guides the user through the testing process step by step to ensure correct usage.

Content

- 1 x 10 parameter dipstick
- 1 x urine receptacle and lid
- 1 x Colour-Board
- User manual
- Smartphone application can easily be downloaded from the app store

